



clinigo[®]

USER MANUAL

Clinical Disinfectant Universal Wipes

Disposable Cleaning Wipes for Disinfecting Surfaces and Non-invasive Medical Devices



Clinical Disinfectant Universal Wipes

50
WIPES

100
WIPES

200
WIPES

IDEAL FOR USE IN HOSPITALS AND HEALTHCARE SETTINGS.

A ready-to-use, single-use, 2-in-1 wipe for the one-step cleaning and disinfection of non-porous surfaces and non-invasive medical devices in a healthcare environment. The wipes are intended for use by healthcare professionals to reduce microbial contamination (including bacteria, yeast and viruses), thereby supporting infection control and prevention protocols.

THE ESSENTIALS

- Bactericidal
- Tuberculocidal
- Mycobactericidal
- Virucidal
- Yeastocidal
- Alcohol free
- Ready to use
- Gentle on surfaces

COMPATIBLE MATERIALS

Compatibility with common healthcare equipment materials such as:

- Polycarbonate
- Acrylonitrile Butadine Styrene (ABS)
- Polymethyl Methacrylate (PMMA) (Acrylic)
- High Density Polyethylene (HDPE)
- Polypropylene
- Stainless Steel

IDEAL FOR

- Effective cleaning and disinfection in a short contact time
- Intermediate disinfection where there is a known or a risk of infection*
- All healthcare environments; available in flow wraps
- Disinfecting surfaces of non-invasive medical devices and equipment

***Not intended for end point sterilisation of medical devices.**

INGREDIENTS

< 5% Non-ionic Surfactants,
< 5% Cationic Surfactants,
Disinfectant, Parfum.

ACTIVE SUBSTANCE PER 100GRAMS OF PRODUCT:


Didecyldimethylammonium chloride (DDAC) 0.86g/100g,
Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-16)) 0.8g/100g,
N-(3-aminopropyl)-N-dodecylpropane-1, 3-diamine (Diamine) 0.7g/100g).

CLINICALLY PROVEN

Proven to inactivate ≥99.99% of tested microorganisms.

*Refer to the efficacy spectrum table for validated organisms and contact times.

Note: All microorganisms listed within the efficacy spectrum table have been tested and demonstrated to achieve at least ≥99.99% (≥4 log reduction) inactivation under the validated test conditions.



Clinical Disinfectant Universal Wipes

50
WIPES

100
WIPES

200
WIPES

DIRECTIONS FOR USE AND CONTRA-INDICATIONS

Use gloves and appropriate PPE. Dispense wipe and open flat to clean the surface, wipe the surface in an 'S' shape motion from top to bottom. To ensure the surface remains wet for the required contact time*, if necessary, use additional wipes. Allow the surface to air dry naturally before use.

CONTRA-INDICATIONS

Do not use wipes for end point sterilisation of medical devices.


STEP ONE

- Use gloves and appropriate PPE
- Follow your organisation's risk assessment policy when using PPE



STEP TWO

- Open lid and label
- Pull out wipe
- Close lid re-sealing pack



STEP THREE

- Clean surface in 'S' shape motion, working from clean to dirty areas
- Keep the wipe flat (not scrunched) to ensure full surface coverage
- Avoid wiping the same area twice



STEP FOUR

- Dispose of wipe in appropriate bin
- Do not re-use wipes
- Allow the surface to dry naturally before use
- Leave the surface for the stated contact time



USE BIOCIDES SAFELY

Always read the label and product information before use.

Re-use of wipes increases the risk of infection, discard each wipe after single use.

FOR PROFESSIONAL USE ONLY

WARNINGS

- Harmful to aquatic life with long lasting effects.
- Keep out of reach of children.
- Avoid release to the environment.
- Dispose of contents/container to a licensed facility in accordance with national regulations.
- Always read the label before use. Use biocides safely and sustainably, this product should be used in accordance with the product label
- Wash hands thoroughly after use.
- If irritation occurs, stop use immediately and rinse with clean water, if symptoms persist seek medical advice
- Do not use with other cleaning products
- Not for skin cleansing or other personal use
- Do not flush down the toilet.
- Do not use the wipes if the packaging is damaged, improperly sealed, leaking, or shows signs of contamination.
- Do not use wipes that appear dry or insufficiently saturated.
- If any visible change in product appearance, odour, colour, or wipe saturation is observed, discontinue use and replace the product.

TO REPORT A POISONING INCIDENT

If a serious incident occurs in relation to this device, the user must report it to the manufacturer and the competent authority of the Member State where they are established.

National Poisons Information Centre
(01 8092166) (IRE)
8am -10pm every day.

For Safety Data Information
(medical personnel only)
T: 08708033979 (UK).

*To achieve the full claimed efficacy spectrum against bacteria, enveloped viruses, yeast and mycobacteria, ensure the surface remains visibly wet for a minimum contact time of 2 minutes (120 seconds).

Clinical Disinfectant Universal Wipes

50
WIPES

100
WIPES

200
WIPES

EFFICACY DETAILS

The liquid extracted from the wipes have been tested against: EN16615, EN14476, EN13727, EN13624 and EN14348 under dirty conditions.

EFFICACY SPECTRUM	TESTING NORM	CONTACT TIME
BACTERICIDAL		
Enterococcus hirae	EN 13727 + EN 16615	30 seconds
Pseudomonas aeruginosa	EN 13727 + EN 16615	30 seconds
Staphylococcus aureus	EN 13727 + EN 16615	30 seconds
Escherichia coli	EN 13727	30 seconds
Klebsiella pneumoniae	EN 13727 + EN 16615	30 seconds
Vancomycin-Resistant Enterococci (VRE)	EN 13727 + EN 16615	30 seconds
Methicillin-resistant Staphylococcus aureus (MRSA)	EN 13727 + EN 16615	30 seconds
YEASTICIDAL		
Candida albicans	EN 13624 + EN 16615	60 seconds
TUBERCULOCIDAL		
Mycobacterium terrae	EN14348	60 seconds
MYCOBACTERICIDAL		
Mycobacterium terrae	EN 14348	60 seconds
Mycobacterium avium	EN 14348	120 seconds
VIRUCIDAL		
Vaccinia virus	EN 14476	15 seconds
Murine Norovirus	EN 14476	30 seconds

*TESTED UNDER DIRTY CONDITIONS




Proven to inactivate ≥99.99% of tested microorganisms.

*Refer to the efficacy spectrum table for validated organisms and contact times.

Note: All microorganisms listed within the efficacy spectrum table have been tested and demonstrated to achieve at least ≥99.99% (≥4 log reduction) inactivation under the validated test conditions.

To achieve the full claimed efficacy spectrum against bacteria, enveloped viruses, yeast, and mycobacteria, ensure the surface remains visibly wet for a minimum contact time of 2 minutes (120 seconds).

PACKAGING SYMBOLS



FOR INFORMATION

For Safety Data Information

(medical personnel only)

T: 08708033979 (UK).



Saraco Industries Ltd

Parkside House
167 Chorley New Road
Bolton, BL1 4RA, United Kingdom
T : +44 (0) 870 803 3979 UK

Clinical Disinfectant Universal Wipes

50
WIPES

100
WIPES

200
WIPES

PRODUCTS WITHIN THIS RANGE

50 Clinical Disinfectant Universal Wipes

Wipe Size:	200 x 200 mm
Substrate:	30gsm
Packs Per Case:	24
Product Code:	CU8676
Barcode:	5060110228676



100 Clinical Disinfectant Universal Wipes

Wipe Size:	200 x 200 mm
Substrate:	30gsm
Packs Per Case:	12
Product Code:	CU8683
Barcode:	5060110228683



200 Clinical Disinfectant Universal Wipes

Wipe Size:	200 x 270 mm
Substrate:	30gsm
Packs Per Case:	6
Product Code:	CU8690
Barcode:	5060110228690



Saraco Industries Ltd
 Parkside House
 167 Chorley New Road
 Bolton, BL1 4RA, United Kingdom
 T : +44 (0) 870 803 3979 UK